

20240812 EB_MN WT-11 DPR, hypo-glu-rescue

Project: Yi-Shuan Tseng Holt Lab
Author: Yi Shuan Tseng
Status: In progress
Entry Created On: 05 Sep 2024 22:10:43 UTC
Entry Last Modified: 05 Sep 2024 22:20:08 UTC
Export Generated On: 01 Sep 2026 18:31:29 UTC

- 2024年8月12日 星期一
- split iPSCs into a full 6-well plate on the Friday before differentiation
differentiate 5 of 6-wells with 10mL T2 in cell-repellent plates
keep some iPSCs
- 2024年8月14日 星期三
- T2 -> T3
- 2024年8月16日 星期五
- T3 -> T4 (~14 mL for over-weekend)
- 2024年8月19日 星期一
- T4 -> T5
- 2024年8月21日 星期三
- T5 -> T6
- 2024年8月23日 星期五
- T6 -> T6 (~14mL for over-weekend)
- 2024年8月26日 星期一
- T6 -> T7
Dissociation
 1. Coat 0.01% Poly-L-Ornithine (PLO) Solution with sodium borate (1:100) for 4-5 hrs
 2. Remove the PLO and wash the plate with PBS once
 3. Coat the plate with laminine for 1hr.
 4. Collect EBs in 50mL falcon and spin them down with 200rcf, 1min. Remove the supernatent (can leave some liquid to avoid remove EBs).
 5. Add trypsin (usually 3-4 mL/10-cm plate) and pipette it up and down with 5mL serological pipette.
 6. Put the tube into 37C water bath. Pipette the cells with P1000 every 5 min round 8-10 times until 15 min. (3 times in total. There might be more resistance at the 2nd round)

7. Add 2x volumn of Heat Inactivated FBS into the tube and pipette it with 5mL serological pipette gently. (Try to make the viscous cluster get through the tip to break it down)
8. Add 2x volumn of CTWM and pipette it with 5mL serological pipette gently.
9. Filter the cell with 40um filter into a new 50mL tube
10. Spin the cell down at 300rcf, 7min. Remove the supernatent.
11. Add T7 (~5mL) to resuspend the cells and count it.
12. Plate the cells with T7. (350-500K/well for 12well plate; 1-1.5M/well for 6-well plate)

Table3								
	A	B	C	D	E	F	G	H
1	cell	passage	trypsin	HIFBS/CTWM	cell count (M/mL)	volume	total	per-12well (K)
2	WT11 (NYU92)	P13	3.5	7	1.16	5	5.8	483

Well1				
	1	2	3	4
A	110 (~>20uL)		111 (~>20uL)	
B				
C				

- 2024年8月27日 星期二
- Change T7 -> T7
- 2024年8月29日 星期四
- T7 -> BP T8
- 2024年8月30日 星期五
- 400uL BP T8 -> 500 uL BP T8
- 2024年9月2日 星期一
- 400uL BP T8 -> 500 uL BP T8
- 2024年9月4日 星期三
- 400uL BP T8 -> 500 uL BP T8
- 2024年9月5日 星期四
- treat cells for condition 6

conditions		^
	A	
1	BP T8 (NT)	
2	1uM PR30	
3	1uM GR30	
4	1uM GP30	
5	hypo 1h	
6	hypo 1h + 24h rescue	

conditions on the plate					^
	1	2	3	4	
A	2	1	4	3	
B	3	4	1	2	
C	5	6	5	6	

treat C2 and C4 with Hypo-Glu (0.15mM glucose) media from 5:30 – 6:30 pm and change it back to BP T8 at 6:30 pm on 9/5

2024年9月6日 星期五

mark positions:
ROWA & ROWB: