

20240826 EB_MN (GEM, TDP43, viability, ATP/ADP)

Project: Yi-Shuan Tseng Holt Lab
Author: Yi Shuan Tseng
Status: In progress
Entry Created On: 16 Sep 2024 18:08:15 UTC
Entry Last Modified: 24 Mar 2025 05:59:27 UTC
Export Generated On: 01 Sep 2026 18:29:42 UTC

- 2024年8月26日 星期一
- split iPSCs into a full 6-well plate on the Friday before differentiation
differentiate 5 of 6-wells with 10mL T2 in cell-repellent plates
keep some iPSCs
- 2024年8月28日 星期三
- T2 -> T3
- 2024年8月30日 星期五
- T3 -> T4 (~14 mL for over-weekend)
- 2024年9月2日 星期一
- T4 -> T5
- 2024年9月4日 星期三
- T5 -> T6
- 2024年9月6日 星期五
- T6 -> T6 (~14mL for over-weekend)
- 2024年9月9日 星期一
- T6 -> T7
Dissociation
 1. Coat 0.01% Poly-L-Ornithine (PLO) Solution with sodium borate (1:100) for 4-5 hrs
 2. Remove the PLO and wash the plate with PBS once
 3. Coat the plate with laminine for 1hr.
 4. Collect EBs in 50mL falcon and spin them down with 200rcf, 1min. Remove the supernatent (can leave some liquid to avoid remove EBs).
 5. Add trypsin (usually 3-4 mL/10-cm plate) and pipette it up and down with 5mL serological pipette.
 6. Put the tube into 37C water bath. Pipette the cells with P1000 every 5 min round 8-10 times until 15 min. (3 times in total. There might be more resistance at the 2nd round)

7. Add 2x volumn of Heat Inactivated FBS into the tube and pipette it with 5mL serological pipette gently. (Try to make the viscous cluster get through the tip to break it down)
8. Add 2x volumn of CTWM and pipette it with 5mL serological pipette gently.
9. Filter the cell with 40um filter into a new 50mL tube
10. Spin the cell down at 300rcf, 7min. Remove the supernatent.
11. Add T7 (~5mL) to resuspend the cells and count it.
12. Plate the cells with T7. (350-500K/well for 12well plate; 1-1.5M/well for 6-well plate)

Table3

	A	B	C	D	E	F	G	H
1	cell	passage	trypsin	HIFBS/CTWM	cell count (M/mL)	volume	total	per 24-well (K)
2	WT11 (NYU92)	P14	3.5+0.5(1st)	7	1.18	5	5.9	246

Well1

	1	2	3	4	5	6
A	51 (cytoGEM)			IF		
B						
C				viability		
D						
	92 (ATP/ADP, TDP43)					

2024年9月10日 星期二

Change T7 -> T7 if added virus

2024年9月11日 星期三

T7-> BP T8 if not added virus

2024年9月12日 星期四

T7 -> BP T8 if added virus

2024年9月13日 星期五

400uL BP T8 -> 500 uL BP T8

2024年9月16日 星期一

400uL BP T8 -> 500 uL BP T8

GEM

	A	B	C	D
1		1. stress	2. 1h resue	3. 24h rescue
2	A. BP			
3	B. hypo 1h			
4	C. hypo 24h			

ATP/ADP, TDP43 (rescue in the same wells)

	A	B	C
1	BP	hypo 1h	hypo 24

IF

	A	B	C	D
1		BP	hypo 1h	hypo24h
2	stess			
3	24h rescue			

viability

	A	B	C	D
1		BP	hypo 1h	hypo24h
2	stess			
3	24h rescue			

2024年9月18日 星期三

400uL BP T8 -> 500 uL BP T8