

20241028EB_MN (qPCR)

Project: Yi-Shuan Tseng Holt Lab
Author: Yi Shuan Tseng
Status: In progress
Entry Created On: 12 Nov 2024 22:38:38 UTC
Entry Last Modified: 24 Mar 2025 06:01:56 UTC
Export Generated On: 01 Sep 2026 18:28:26 UTC

- 2024年10月28日 星期一
- split iPSCs into a full 6-well plate on the Friday before differentiation
differentiate 5 of 6-wells with 10mL T2 in cell-repellent plates
keep some iPSCs
- 2024年10月30日 星期三
- T2 -> T3
- 2024年11月1日 星期五
- T3 -> T4 (~14 mL for over-weekend)
- 2024年11月4日 星期一
- T4 -> T5
- 2024年11月6日 星期三
- T5 -> T6
- 2024年11月8日 星期五
- T6 -> T6 (~14mL for over-weekend)
- 2024年11月11日 星期一
- T6 -> T7
- Dissociation
1. Coat 0.01% Poly-L-Ornithine (PLO) Solution with sodium borate (1:100) for 4-5 hrs
 2. Remove the PLO and wash the plate with PBS once
 3. Coat the plate with laminine for 1hr.
 4. Collect EBs in 50mL falcon and spin them down with 200rcf, 1min. Remove the supernatent (can leave some liquid to avoid remove EBs).
 5. Add trypsin (usually 3-4 mL/10-cm plate) and pipette it up and down with 5mL serological pipette.
 6. Put the tube into 37C water bath. Pipette the cells with P1000 every 5 min round 8-10 times until 15 min. (3 times in total. There might be more resistance at the 2nd round)
 7. Add 2x volumn of Heat Inactivated FBS into the tube and pipette it with 5mL serological pipette gently. (Try to make the viscous cluster get through the tip to break it down)
 8. Add 2x volumn of CTWM and pipette it with 5mL serological pipette gently.

9. Filter the cell with 40um filter into a new 50mL tube
10. Spin the cell down at 300rcf, 7min. Remove the supernatent.
11. Add T7 (~5mL) to resuspend the cells and count it.
12. Plate the cells with T7. (350-500K/well for 12well plate; 1-1.5M/well for 6-well plate)

put 2 plates into 2 tubes (did not work v. well...)

Table1 ^

	A	B
1		mL/tube
2	trypsin	3.5+0.5
3	HI-FBS/CTWM	8

Table2 ^

	A	B	C	D	E
1	conc. (M/mL)	survival rate	total volume	total cell #	cell#/well
2	1.44	97%	5	7.2	1.2

plate cells into 1 full 6-well dish for QPCR

Well1 ^

	1	2	3
A	2.5mM 1h	0.15mM 1h	0.15mM 24h
B	2.5mM 1h + 24h rescue	0.15mM 1h + 24h rescue	0.15mM 24h + 24h rescue

2024年11月13日 星期三

T7-> BP T8

2024年11月15日 星期五

400uL BP T8 -> 500 uL BP T8

2024年11月18日 星期一

400uL BP T8 -> 500 uL BP T8

2024年11月20日 星期三

400uL BP T8 -> 500 uL BP T8 for A1-3 and B1-2

treat B3 with hypo at aou