2025 iPSC record

Project: Yi-Shuan Tseng Holt Lab
Author: Yi Shuan Tseng
Status: In progress
Entry Created On: 17 Jun 2025 21:40:27 UTC
Entry Last Modified: 18 May 2026 19:15:26 UTC
Export Generated On: 01 Sep 2026 18:25:26 UTC

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Thaw new iPSC from liquid nitrogen
hiPSC NIH wt-11 P16 Halo-Tag-TARDBP (N-term) 1/2 6-well 20240926 YST
into 1* 6-well in STEMFlex + ROCKI
- 2025年6月18日 星期三

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The plate was very full (too confluent)
Differentiate 5 wells and split 1 well to 2 well (1:6) (P2)
I accidentally used BP T8 to split the cells --> collect the cells --> spin it down --> wash it with Stemflex --> spin down --> resuspend in stemflex
- 2025年6月25日 星期三

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Found cells have yeast (maybe pombe) contamination
Toss the cells and investigate media and matrigel

Well1			
	1	2	3
A	matrigel + Stemflex+ROCK I	matrigel+ Stemflex	matrigel
B	Stemflex+ROCK I	Stemflex	

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A1 looked a bit contaminated
Thaw new cells with new matrigel and aliquot media from stock. I used new Corning DMEM to prepare it.

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Change STEMFLEX + ROCKI to STEMFLEX

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Chagne media

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Split cells 1:8

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Change media (4mL) for over-weekend

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Cells are confluent
split 1:8 again (grow much faster than before)
differentiate neurons: 20250818MN (the yield was bad, EBs were very hard to break)

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Passage 1:6

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Cells look normal

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Cells are very confluent and many cells looked differentiated or bigger than usual iPSC
Split 1:8 into 6*6-well

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Chao helped me change media for iPSCs

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Cells looked even more different and super confluent --> trash
Thaw a new vial (hiPSC NIH wt-11 P16 Halo-Tag-TARDBP (N-term) 1/2 6-well 20240926 YST into 2*6-well)

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Split iPSCs into 2* full 6-well (1:8)

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Differentiate 5 wells of 1 plate into T2 [20250908MN](#) (very hard to break EBs. Very low yield)

Passage 1 well into 2 (1:8)

Froze the other plate into 12 vials (

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Cells is almost confluent. Split cells into 6*6-well (1:8) (Plate A)

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Cells grew very fast. I suspected that it might related to the new DMEM for Matrigel since cells grew much faster after I used Corning DMEM to prepare Matrigel.

Use Gibco DMEM +glutaMAX to prepare a new Matrigel coating solution. Thaw a new vial into 2*6-well. (hiPSC NIH wt-11 P16 Halo-Tag-TARDBP (N-term) 1/2 6-well 20240926 YST into 2*6-well) PlateB

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Cells on Plate A are almost confluent. Split one wells into 6*6-well coated Matrigel with Gibco DMEM (1:8) (Plate C)

Plate B: change media

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Plate A is super confluent

Plate B is around 70% confluent -> split 1:6 into 2 6-well plates

Plate C grew slow again. It's only around 30% confluent, but the cell morphology looked much better and close to what I had before.

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Differentiate all 6 wells of Plate C today. [20250916MN](#) (very very hard to break EBs. The lowest yield since I started to use accutase method)

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The cells looked much normal than before. I supposed the abnormal phenotype is related to the DMEM with Corning media.

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Try to use 'tapping' way to detach iPSCs after RelesR according to the protocol.

1. add 1mL of RelesR
2. Aspirate out the RelesR
3. incubate in 37C for 5-7min (may vary for different iPSCs)
4. add 1mL media to the well
5. Tapp the well for 30-60s to detach the undifferentiated iPSCs

The first attempt I incubate the cells for around 4-5min and the most of them detached.

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The cells were around 50% confluent. Split the cells with tapping method of ReleSR. 1:12. (Plate1)

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Split a 70-80% confluent cells 1:14 *3 and 1:18*3.

Plate2				^
	1	2	3	
A	1:14			
B	1:18			

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Plate1 is a round 70% confluent, but there are very big colony. The cells in the middle looked not healthy.
Plate2 is around 70-80% confluent. 1:14 ones are a bit more. --> differentiate 5 wells.
Split B2 of Plate2 into 1:20 (6*6-well --> Plate3) using tapping method of ReleSR. I incubated for 8-9min. Around 60-70% cells were detached. I pipetted the cells with 5mL pipette maybe around 10 times to reach to smaller size clumps. After plating, the size of the clumps are ok (smaller than what I had before, which I thought was a bit too big).

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Plate3 is around 5%

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Since the EB dissociation kept having problem --> Thaw a new vial of WT-11 into 2*6-well
hiPSC NIH wt-11 P16 Halo-Tag-TARDBP (N-term) 1/2 6-well 20250721 YST0003

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Change media

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70-80% confluent
Split all cells into 1 full 6-well plate (1:3) with Stemflex with ROCKI

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Change media to Stemflex

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~60-70% confluent
--> start differentiation
--> keep 1 well (1:3) using accutase