2026 iPSC record (Jan-Jun)

Project: Yi-Shuan Tseng Holt Lab

Author: Yi Shuan Tseng

Status: In progress

Entry Created On: 20 Jan 2026 16:37:19 UTC

Entry Last Modified: 26 Aug 2026 21:10:03 UTC

Export Generated On: 01 Sep 2026 18:34:21 UTC

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Coat 6-well plastic plate with Matrigel (1:100 in DMEM, prepared today) for 1h in 37C

Thaw a new vial of WT-11 into 2*6-well with Stemflex+ROCKI

hiPSC NIH wt-11 P16 Halo-Tag-TARDBP (N-term) 1/2 6-well 20240926 YST

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Cells are around 40-50% confluent

Change media to Stemflex

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Cells is around 90% confluent

Split cells 1:8 into 8 wells

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Cells is around 90% confluent

Split cells 1:8 into 2 6-wells

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Cells is around 70% confluent. There is a big cluster in the middle. May need to break it down more.

Split cells 1:8 into 2 6-wells

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Cells is around 50% confluent.

Split cells with ReleSR. Pipetted more to break the clusters.

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Cells look a bit differentiated (cells was P6). Change media and check them again tomorrow.

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Cells look differentiated. I might pipette them too much? Do not use them next week.

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Thaw cells into 2*6-wells with Stemflex+ROCKI

hiPSC NIH wt-11 P16 Halo-Tag-TARDBP (N-term) 1/2 6-well 20250721 YST0002

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~30% confluent, change media to Stemflex

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Change media. ~50%

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~80%
Split cells with accutase. 180K/well

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Change media to Stemflex

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~60% confluent

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~95% confluent
Split with accutase today with Stemflex + RockI
conc. 494K/mL
200K: 445uL
250K: 506uL

Well1			
	1	2	3
A			200K
B			250K

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The cells didn't grow very well.
The confluency is less than last time when I plated 180K or 200K. (only a half)
Change media

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Cells ~40-50% confluent

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Cells grew ~ 90% -> split 1:5 into 10 wells of 6-well with ReleSR. I didn't pipette too much to avoid over-pipetting. Clusters are big. I tried to shake it evenly.

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Cells grew a lot and confluent in the middle of the well of a plate I coated yesterday (newly diluted matrigel). The rest of the well is very sparse and empty. The other four wells grew relatively evenly (the plate I coated around Thursday or Friday), but it's around 50% confluent. --> differentiate tomorrow.

However, I saw some differentiated-like cells again (It was only P4). Not all the cells, but I saw a bit more then just around the edge. I am pretty confused. I guess maybe the Stemflex supplement was expired too long (I remembered it expired in March 2024). Need to observe again tomorrow.

Therefore, I thaw new cells with freshly thawed ROCKI into 2*6-well.

hiPSC NIH wt-11 P17 Halo-Tag-TARDBP (N-term) 1/2 6-well 20250911 YST0004

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Split cells with accutase today (1:6)

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Change media to Stemflex

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Split iPSC with ReleSR. Only a few cells detached...

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Cells did not survive well. Need to thaw new cells again.

hiPSC NIH wt-11 P17 Halo-Tag-TARDBP (N-term) 1/2 6-well 20250911 YST0005

Thaw

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Did the 5th passage today.

All the passages looked good with newly ordered StemFlex media. I might just use the newly-ordered media from now on.

I didn't have evidence to show the old media has problem though, but cells usually started looking weird by the 4th passage.

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(I dumped the old media at the end since the cells in new media seemed to be much more stable)

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Thaw hiPSC NIH wt-11 P17 Halo-Tag-TARDBP (N-term) 1/2 6-well 20250911 YST0006

into 2 wells of 6-well plate.

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Change into media without ROCKi

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Split cells from 1 well into 10 wells of 6-well with ReleSR

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~10-20% confluency

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30-40% confluency, change media

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~70% confluency, more confluent in the middle and still quite some space around around the edge

Note: I should expand some cells from Martina's storage again to get some earlier passage cells.

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Split cells with ReleSR. Not many cells were detached. I might pipette with P1000 a bit too much. Clusters were broken down into small pieces.

Seed 1:6 and 1:8 into 2 wells of 6-well

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Cells didn't attached well. Need time to grow

At the same time, I thaw a WT11 cells from Martina's Green LN Box into 2 wells of 6-well with STEMFlex + ROCK I (I feel the cells I thawed was not easy to detach by ReLeSR might indicate a more differentiated state.)

2023/1/20 hiPSCs wt-11 p17 N-term HaloTag TARDBP

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Thaw hiPSC NIH wt-11 P19 Halo-Tag-TARDBP (N-term) 1/2 6-well 20260526 YST0124