

# 20260525EB\_MN (contaminated)

**Project:** Yi-Shuan Tseng Holt Lab

**Author:** Yi Shuan Tseng

**Status:** In progress

**Entry Created On:** 11 Jun 2026 15:38:50 UTC

**Entry Last Modified:** 11 Jun 2026 15:48:16 UTC

**Export Generated On:** 01 Sep 2026 18:34:47 UTC

**2026年5月26日 星期二**

split iPSCs into a full 6-well plate on the Friday before differentiation

differentiate 5 of 6-wells with 10mL T2 in cell-repellent plates

keep some iPSCs

**2026年5月28日 星期四**

T2 -> T3

**2026年5月30日 星期六**

T3 -> T4 (~14 mL for over-weekend)

**2026年6月2日 星期二**

T4 -> T5

**2026年6月4日 星期四**

T5 -> T6

**2026年6月6日 星期六**

T6 -> T6 (~14mL for over-weekend)

**2026年6月8日 星期一**

T6 -> T7

Dissociation

1. Coat 0.01% Poly-L-Ornithine (PLO) Solution with sodium borate (1:100) for 4-5 hrs
2. Remove the PLO and wash the plate with PBS once
3. Coat the plate with laminine for 1hr.
4. Collect EBs in 50mL falcon and spin them down with 200rcf, 1min. Remove the supernatent (can leave some liquid to avoid remove EBs).
5. Add trypsin (usually 3-4 mL/10-cm plate) and pipette it up and down with 5mL serological pipette.
6. Put the tube into 37C water bath. Pipette the cells with P1000 every 5 min round 8-10 times until 15 min. (3 times in total. There might be more resistance at the 2nd round)
7. Add 2x volumn of Heat Inactivated FBS into the tube and pipette it with 5mL serological pipette gently. (Try to make the viscous cluster get through the tip to break it down)
8. Add 2x volumn of CTWM and pipette it with 5mL serological pipette gently.

9. Filter the cell with 40um filter into a new 50mL tube
10. Spin the cell down at 300rcf, 7min. Remove the supernatent.
11. Add T7 (~5mL) to resuspend the cells and count it.
12. Plate the cells with T7. (350-500K/well for 12well plate; 1-1.5M/well for 6-well plate)

| cell info |      |         | ^ |
|-----------|------|---------|---|
|           | A    | B       |   |
| 1         | cell | passage |   |
| 2         | WT11 |         |   |
|           |      |         |   |

| Table1 |             |                           | ^ |
|--------|-------------|---------------------------|---|
|        | A           | B                         |   |
| 1      | run1        | volume (mL)               |   |
| 2      | trypsin     | 3                         |   |
| 3      | HI-FBS/CTWM | 6                         |   |
| 4      | run2        | *take some undissolved EB |   |
| 5      | trypsin     | 2                         |   |
| 6      | HI-FBS/CTWM | 4                         |   |

| Table3 |              |          |              |                 |        |            |                     |                 |          | / |
|--------|--------------|----------|--------------|-----------------|--------|------------|---------------------|-----------------|----------|---|
|        | A            | B        | C            | D               | E      | F          | G                   | H               | I        |   |
| 1      | conc. (M/mL) | survival | total volume | tota cell count | plate# | well plate | cell count/well (K) | number of wells | volume   |   |
| 2      | 4.26         | 99       | 7            | 29.82           | A      | 96         | 60                  | 36              | 507      |   |
| 3      |              |          |              |                 | B      | 24         | 250                 | 16              | 939      |   |
| 4      |              |          |              |                 | C,D    | 6          | 1972                | 12              | the rest |   |

2026年6月9日 星期二

Plate I & II were contaminated

2026年6月10日 星期三

Plate III and IV were contaminated

as well as newly started 20260608EB\_MN and PLO coated spheres in Yoshiho's device (Matrigel coated ones were clear)

Suspected that the PBS got contaminated since PBS is the only common thing for all of them

Culture the [PBS, T2, T3, T5, T6, BPT8] with U2OS cells, in empty plate, and in LB borth (I took a colony from a random plate in cold room (perhaps from Andy for a positive control)

2026年6月11日 星期四

Somehow the potential positive control did not grow, but PBS definitely grew.

 PXL\_20260611\_153348786.jpg



The plate with cells can see floating bacteria with PBS, but they haven't exploded

The media only plate. Clear, Opaque with contamination

| Well1 |                         |      |             |   |   |   | ^ |
|-------|-------------------------|------|-------------|---|---|---|---|
|       | 1                       | 2    | 3           | 4 | 5 | 6 |   |
| A     | T6                      | T2   |             |   |   |   |   |
| B     | T6<br>+vLH<br>92<br>1uL | T3   |             |   |   |   |   |
| C     | T6<br>+vLH<br>92<br>5uL | T5   |             |   |   |   |   |
| D     | PBS                     | BPT8 | PBS<br>+ T6 |   |   |   |   |