

20260602EB_MN (NCT, ATP, IF)

Project: Yi-Shuan Tseng Holt Lab

Author: Yi Shuan Tseng

Status: In progress

Entry Created On: 24 Jun 2026 01:40:58 UTC

Entry Last Modified: 08 Jul 2026 22:02:57 UTC

Export Generated On: 01 Sep 2026 18:36:33 UTC

2026年6月2日 星期二

split iPSCs into a full 6-well plate on the Friday before differentiation

differentiate 5 of 6-wells with 10mL T2 in cell-repellent plates

keep some iPSCs

2026年6月4日 星期四

T2 -> T3

2026年6月6日 星期六

T3 -> T4

2026年6月8日 星期一

+3mL T4

2026年6月9日 星期二

T4 -> T5

2026年6月11日 星期四

T5 ->

2026年6月13日 星期六

T6 -> T6 (~14mL for over-weekend)

2026年6月16日 星期二

T6 -> T7

Dissociation

1. Collect EBs into 50mL tubes (each plate). Wash the plate with 5mL PBS and added them into the tubes. Let EBs settle with gravity.
2. Remove sup. and add 3mL 0.25% trypsin for each (prepare 8mL 0.25% trypsin with 220U DNaseI)
3. Incubate the tubes in 37C.
4. pipette*10 at 0min, 5min, 10min.
5. Stop reaction at ~15 min by adding 8mL FBS into 50mL tube -> pipette with 5mL pipet ~ 10 times
6. add 8mL CTWM -> pipette ~ 10 times
7. filter cells through 40um cell strainer into new tubes.
8. centrifuge at 200xg, 4min

cell info ^

	A	B	C	D
1		cell	passage	well#
2	A	WT11	P17+4	3
3	B	WT11	P17+4	4

Table1 ^

	A	B
1		volume (mL)
2	trypsin	4
3	HI-FBS/CTWM	8

Table3 ^

	A	B	C	D	E	F	G	H	I	J
1		conc. (M/mL)	survival	total volume	tota cell count	plate#	well plate	cell count/well (K)	number of wells	volume
2	A	3.09	99	3	9.27	I (A)	24	300	5	485
3	B	3.44	100	3	10.32	I (B)	24	300	5	436
4						II (A) 92	24	300	10	971
5						II (B) 92	24	300	10	872
6						I (A) IF	24	300	2	194
7						I (B) IF	24	300	2	174
8						II (B) 51	24	300	4	349

Batch A

Batch B

Plate I ^

	1	2	3	4	5	6
A	110 (7.5uL, 1.5uL/well)					no virus
B	111 (8uL, 1.6uL/well)					
C	110 (7.5uL, 1.5uL/well)					no virus
D	111 (8uL, 1.6uL/well)					

Plate II						
	1	2	3	4	5	6
A						51 1uL
B	92 (12,5uL for 10wells)					51 1uL
C						51 2uL
D	92 (12,5uL for 10wells)					51 2uL

- 2026年6月17日 星期三
- T7-> T7
- 2026年6月19日 星期五
- T7 -> T8
- 2026年6月21日 星期日
- 1/2 T8 -> 1/2 T8
- 2026年6月23日 星期二
- 1/2 T8 -> 1/2 T8 (new media (unexpired))
- 2026年6月24日 星期三

Imaging vLH92 (batch A on plate II) - MMS treatment

- ** Cells looked better than last batch
- Do primary dilution and then do the secondary dilution right before imaging. Imaging from E -> A
1. Prepare 50*1 2.5mM gluc NBA + T8 supp + 1mM pyruvate

2. Aliquot 20uL MMS into eppendorf and cover it with tin foil in RT

3. Dilute MMS 1uL in 99uL media right before the treatment

4. Prepare different concentration of MMS

5. Mix well

6. Right before doing olaparib experiment

a. Dilute olaparib (10mM stock) 1:10 to 1mM (1uL into 9uL media)

b. Take 11mL media and add 2.2uL 1mM olaparib

Plate II condition						
	1	2	3	4	5	6
A	A	B	C	D	E	
B	A'	B'	C'	D'	E'	
C	D	C	B	A	E	
D	D'	C'	B'	A'	E'	

condition1								
	A	B	C	D	E	F	G	H
1	condition #	MMS	media		condition #	MMS	treatment	media
2	A	MMS 1.5mM	with 1mM pyruvate		A'	MMS 1.5mM	olaparib 200nM	with 1mM pyruvate
3	B	MMS 1mM	with 1mM pyruvate		B'	MMS 1mM	olaparib 200nM	with 1mM pyruvate
4	C	MMS 0.5mM	with 1mM pyruvate		C'	MMS 0.5mM	olaparib 200nM	with 1mM pyruvate
5	D	MMS 0.25mM	with 1mM pyruvate		D'	MMS 0.25mM	olaparib 200nM	with 1mM pyruvate
6	E	MMS 0mM	with 1mM pyruvate		E'	MMS 0mM	olaparib 200nM	with 1mM pyruvate

prepare media with different conc. of glucose with pyruvate1							
	A	B	C	D	E	F	G
1	diluted from (mM)	targetconc. (mM)	dilution factor	final volume (uL)	added volume (uL)	media volume	100mM pyruvate (uL)
2	1100	2.5	440	50000	113.6	49886.4	500
3			total media volume	50.0	mL		

drug calculation													
	A	B	C	D	E	F	G	H	I	J	K	L	M
1			stock (mM)	1st dilution factor	total volume for 1st dil.	stock volume	media volume for 1st dil.	conc. after 1st dilution (mM)	working (mM)	dilution factor from 1st dil.	total volume (mL)	volume of 1st dil.	volume of media
2	MMS	E							0		1000	0	1000
3		D	11800	100	100	1	99	118	0.25	472	1000	2.1	997.9
4		C	11800	100	100	1	99	118	0.5	236	1000	4.2	995.8
5		B	11800	100	100	1	99	118	1	118	1000	8.5	991.5
6		A	11800	100	100	1	99	118	2	59	1000	16.9	983.1
7										total	5000	31.8	

time table for Perceval treatment_MMS plate I											
	A	B	C	D	E	F	G	H	I	J	K
1	Perceval		Plate I								
2	condition#		condition	drug	media	pyruvate	time	imaging	washout with 2.5mM gluc	note	recovery observation on 6/25 with bench microscope in TC
3											similar to pre-imaging
4	E	A5	MMS 0mM	NT	2.5mM gluc	1mM	11:40	12:45			dead
5	A	A1	MMS 2mM	NT	2.5mM gluc	1mM	12:03	13:05			dead
6	B	A2	MMS 1mM	NT	2.5mM gluc	1mM	12:30	13:30			a lot of dead cells
7	C	A3	MMS 0.5mM	NT	2.5mM gluc	1mM	12:50	13:50			some cells survived, but many dead cells
8	D	A4	MMS 0.25mM	NT	2.5mM gluc	1mM	13:10	14:12			
9											similar to pre-imaging
10	E'	B5	MMS 0mM	olaparib 200nM	2.5mM gluc	1mM	13:37	14:35			some of them died, but some were fine
11	A'	B1	MMS 2mM	olaparib 200nM	2.5mM gluc	1mM	13:55	14:55			similar to pre-imaging
12	B'	B2	MMS 1mM	olaparib 200nM	2.5mM gluc	1mM	14:15	15:16			similar to pre-imaging
13	C'	B3	MMS 0.5mM	olaparib 200nM	2.5mM gluc	1mM	14:40	15:40			
14	D'	B4	MMS 0.25mM	olaparib 200nM	2.5mM gluc	1mM	15:00	16:01	right after the next treatment		

fileid_treat_20260624_MN_92						
	A	B	C	D	E	F
1	image	treatment	drug	treat_time	pyruvate	replicate
2	vLH92_MN_60x-1xzoom_0000	0mM MMS	NT	1h	yes	A
3	vLH92_MN_60x-1xzoom_0001	2mM MMS	NT	1h	yes	A
4	vLH92_MN_60x-1xzoom_0002	1mM MMS	NT	1h	yes	A
5	vLH92_MN_60x-1xzoom_0003	0.5mM MMS	NT	1h	yes	A
6	vLH92_MN_60x-1xzoom_0004	0.25mM MMS	NT	1h	yes	A
7	vLH92_MN_60x-1xzoom_0005	0mM MMS	200nM olaparib	1h	yes	A
8	vLH92_MN_60x-1xzoom_0006	2mM MMS	200nM olaparib	1h	yes	A
9	vLH92_MN_60x-1xzoom_0007	1mM MMS	200nM olaparib	1h	yes	A
10	vLH92_MN_60x-1xzoom_0008	0.5mM MMS	200nM olaparib	1h	yes	A
11	vLH92_MN_60x-1xzoom_0009	0.25mM MMS	200nM olaparib	1h	yes	A

2026年6月25日 星期四

Imaging vLH92 (batch B on plate II) - MMS treatment

** Note: cells looked great

Do primary dilution and then do the secondary dilution right before imaging. Imaging from E -> A

1. Prepare 50*1 2.5mM gluc NBA + T8 supp + 1mM pyruvate
2. Aliquot 20uL MMS into eppendorf
3. Dilute MMS 1uL in 99uL media right before the experiment
4. Prepare different concentration of MMS
5. Mix well
6. Right before doing olaparib experiment

a. Dilute olaparib (10mM stock) 1:10 to 1mM (1uL into 9uL media)

b. Take 5mL media and add 1uL 1mM olaparib

Plate II condition1						
	1	2	3	4	5	6
A	A	B	C	D	E	
B	A'	B'	C'	D'	E'	
C	D	C	B	A	E	
D	D'	C'	B'	A'	E'	

condition2								
	A	B	C	D	E	F	G	H
1	condition #	MMS	media		condition #	MMS	treatment	media
2	A	MMS 1.5mM	with 1mM pyruvate		A'	MMS 1.5mM	olaparib 200nM	with 1mM pyruvate
3	B	MMS 1mM	with 1mM pyruvate		B'	MMS 1mM	olaparib 200nM	with 1mM pyruvate
4	C	MMS 0.5mM	with 1mM pyruvate		C'	MMS 0.5mM	olaparib 200nM	with 1mM pyruvate
5	D	MMS 0.25mM	with 1mM pyruvate		D'	MMS 0.25mM	olaparib 200nM	with 1mM pyruvate
6	E	MMS 0mM	with 1mM pyruvate		E'	MMS 0mM	olaparib 200nM	with 1mM pyruvate

prepare media with different conc. of glucose with pyruvate2							
	A	B	C	D	E	F	G
1	diluted from (mM)	targetconc. (mM)	dilution factor	final volume (uL)	added volume (uL)	media volume	100mM pyruvate (uL)
2	1100	2.5	440	50000	113.6	49886.4	500
3			total media volume	50.0	mL		

drug calculation1													
	A	B	C	D	E	F	G	H	I	J	K	L	M
1			stock (mM)	1st dilution factor	total volume for 1st dil.	stock volume	media volume for 1st dil.	conc. after 1st dilution (mM)	working (mM)	dilution factor from 1st dil.	total volume (mL)	volume of 1st dil.	volume of media
2	MM S	E							0		1000	0	1000
3		D	11800	100	100	1	99	118	0.25	472	1000	2.1	997.9
4		C	11800	100	100	1	99	118	0.5	236	2500	10.6	2489.4
5		B	11800	100	100	1	99	118	1	118	2500	21.2	2478.8
6		A	11800	100	100	1	99	118	2	59	2500	42.4	2457.6
7										total	9500	76.3	

time table for Perceval treatment_MMS plate I1															
	A	B	C	D	E	F	G	H	I	J	K	L	M	N	
1	Perceval		Plate I												
2	condition#		condition	drug	media	pyruvate	imaging order	time	imaging	washout with 2.5mM gluc	note	recovery observation on 6/26 with bench microscope in TC			
3	NT														
4	E	C 5	MMS 0mM	NT	2.5mM gluc	1mM	1	17:20	18:25			similar to pre-imaging			
5	A	C 1	MMS 2mM	NT	2.5mM gluc	1mM	7	19:38	20:41			dead			
6	B	C 2	MMS 1mM	NT	2.5mM gluc	1mM	8	20:00	21:02			dead			
7	C	C 3	MMS 0.5mM	NT	2.5mM gluc	1mM	9	20:27	21:27			a lot of dead cells			
8	D	C 4	MMS 0.25mM	NT	2.5mM gluc	1mM	10	20:50	21:50			some cells survived, but many dead cells			
9	200nM olaparib														
10	E'	D 5	MMS 0mM	olaparib 200nM	2.5mM gluc	1mM	2	17:45	20:25			similar to pre-imaging			
11	A'	D 6	MMS 2mM	olaparib 200nM	2.5mM gluc	1mM	3	18:10	19:12			some of them died, but some were fine			
12	B'	D 7	MMS 1mM	olaparib 200nM	2.5mM gluc	1mM	4	18:30	19:30			similar to pre-imaging			
13	C'	D 8	MMS 0.5mM	olaparib 200nM	2.5mM gluc	1mM	5	18:50	19:51			similar to pre-imaging			
14	D'	D 9	MMS 0.25mM	olaparib 200nM	2.5mM gluc	1mM	6	19:17	20:17	right after the next treatment		similar to pre-imaging			

fileid_treat_20260625_MN_92						
	A	B	C	D	E	F
1	image	treatment	drug	treat_time	pyruvate	replicate
2	vLH92_MN_60x-1xzoom_0000	0mM MMS	NT	1h	yes	B
3	vLH92_MN_60x-1xzoom_0001	0mM MMS	200nM olaparib	1h	yes	B
4	vLH92_MN_60x-1xzoom_0002	2mM MMS	200nM olaparib	1h	yes	B
5	vLH92_MN_60x-1xzoom_0003	1mM MMS	200nM olaparib	1h	yes	B
6	vLH92_MN_60x-1xzoom_0004	0.5mM MMS	200nM olaparib	1h	yes	B
7	vLH92_MN_60x-1xzoom_0005	0.25mM MMS	200nM olaparib	1h	yes	B
8	vLH92_MN_60x-1xzoom_0006	2mM MMS	NT	1h	yes	B
9	vLH92_MN_60x-1xzoom_0007	1mM MMS	NT	1h	yes	B
10	vLH92_MN_60x-1xzoom_0008	0.5mM MMS	NT	1h	yes	B
11	vLH92_MN_60x-1xzoom_0009	0.25mM MMS	NT	1h	yes	B

2026年6月26日 星期五

Imaging NCT: Plate I

SirDNA: 1:6000, 1.75uL SirDNA in 10.5mL 2.5mM gluc + 1mM pyruvate
Did not wash out MMS -> just aspirated media after imaging were all finished

NCT_O.C. , 60X, 600x900, 16-bit					
	A	B	C	D	E
1	channel	bit	power (%)	exposure time	note
2	561	16	25 (111)/ 30 (110)	300 ms	
3	640	16	10	100 ms	
4	488	16	50	1s	

time table for NCT treatment_MMS_Plate I-BatchA											
	A	B	C	D	E	F	G	H	I	J	K
1	NCT (plate I)										
2	condition#		condition	drug	media	pyruvate	imaging order	time	imaging	note	recovery observation on 6/14
3		110									
4	E	A5	MMS 0mM	NT	2.5mM gluc	1mM	1	10:15	11:17	cells did not look very good	
5	A	A1	MMS 2mM	NT	2.5mM gluc	1mM	3	11:38	12:36	cells look okay (better than A5)	
6	B	A2	MMS 1mM	NT	2.5mM gluc	1mM	5	12:56	13:55		
7	C	A3	MMS 0.5mM	NT	2.5mM gluc	1mM	7	14:15	15:15	a bit harder to find cells clustering together, but better than A5	
8	D	A4	MMS 0.25mM	NT	2.5mM gluc	1mM	9	17:30	18:28	position 4 shifted focus at the end of the movie	
9		111									
10	E	B5	MMS 0mM	NT	2.5mM gluc	1mM	2	10:55	11:53	cells looked good (better than 110)	
11	A	B1	MMS 2mM	NT	2.5mM gluc	1mM	4	12:15	13:13		
12	B	B2	MMS 1mM	NT	2.5mM gluc	1mM	6	13:36	14:35		
13	C	B3	MMS 0.5mM	NT	2.5mM gluc	1mM	8	16:50	17:50		
14	D	B4	MMS 0.25mM	NT	2.5mM gluc	1mM	10	18:10	19:10		

fileID_treatment_NCT_20260626_MN_110,111						
	A	B	C	D	E	F
1	image	reporter	treatment	drug	time	pyruvate
2	com_110_MN_no-glu_50%_1s_30s_60x_0000	110	0mM MMS	NT	1h	yes
3	com_110_MN_no-glu_50%_1s_30s_60x_0001	110	0mM MMS	NT	1h	yes
4	com_110_MN_no-glu_50%_1s_30s_60x_0002	110	2mM MMS	NT	1h	yes
5	com_110_MN_no-glu_50%_1s_30s_60x_0003	110	2mM MMS	NT	1h	yes
6	com_110_MN_no-glu_50%_1s_30s_60x_0004	110	1mM MMS	NT	1h	yes
7	com_110_MN_no-glu_50%_1s_30s_60x_0005	110	1mM MMS	NT	1h	yes
8	com_110_MN_no-glu_50%_1s_30s_60x_0006	110	0.5mM MMS	NT	1h	yes
9	com_110_MN_no-glu_50%_1s_30s_60x_0007	110	0.5mM MMS	NT	1h	yes
10	com_110_MN_no-glu_50%_1s_30s_60x_0008	110	0.25mM MMS	NT	1h	yes
11	com_110_MN_no-glu_50%_1s_30s_60x_0009	110	0.25mM MMS	NT	1h	yes
12	com_111_MN_no-glu_50%_1s_30s_60x_0000	111	0mM MMS	NT	1h	yes
13	com_111_MN_no-glu_50%_1s_30s_60x_0001	111	0mM MMS	NT	1h	yes
14	com_111_MN_no-glu_50%_1s_30s_60x_0002	111	2mM MMS	NT	1h	yes
15	com_111_MN_no-glu_50%_1s_30s_60x_0003	111	2mM MMS	NT	1h	yes
16	com_111_MN_no-glu_50%_1s_30s_60x_0004	111	1mM MMS	NT	1h	yes
17	com_111_MN_no-glu_50%_1s_30s_60x_0005	111	1mM MMS	NT	1h	yes
18	com_111_MN_no-glu_50%_1s_30s_60x_0006	111	0.5mM MMS	NT	1h	yes
19	com_111_MN_no-glu_50%_1s_30s_60x_0007	111	0.5mM MMS	NT	1h	yes
20	com_111_MN_no-glu_50%_1s_30s_60x_0008	111	0.25mM MMS	NT	1h	yes
21	com_111_MN_no-glu_50%_1s_30s_60x_0009	111	0.25mM MMS	NT	1h	yes

Imaging PfV cytoGEMs (vLH51)

Plate II condition_51						
	1	2	3	4	5	6
A						1
B						4
C						2
D						3

condition					^
	A	B	C	D	
1	condition#	gluc conc (mM)	media		
2	1	2.5	1mM pyruvate		
3	2	0.15	1mM pyruvate		
4	3	0	1mM pyruvate		
5	4	0	1mM pyruvate	+6mM 2DG + 1uM FCCP	

prepare media with different conc. of glucose with pyruvate								^
	A	B	C	D	E	F	G	
1		diluted from (mM)	targetconc. (mM)	dilution factor	final volume (uL)	added volume (uL)	media volume	
2	2	2.5	0.15	16.6666666667	1000	60.0	940.0	
3				total media volume	1.0	mL		
4	2DG	300	6	50	1000	20.0	980.0	
5	FCCP	1	0.001	1000	1000	1.0	999.0	

Treat_time_vLH51													✓
	A	B	C	D	E	F	G	H	I	J	K	L	
1					vLH51								
2	condition#		condition	Halo dye concentration	treat time + Halo_dye	image time	soma	neurite	TDP43				
3	1	A6	2.5mM glucose 1h	10nM	18:40	19:40		1-15	1-10	not totally sure if incubator light is on during TDP43Halo	the GEMs in soma are mostly too much, but can still find cells with good number of GEMs		
4	2	C6	0.15mM glucose 1h	10nM	19:10	20:10		16-25	11-20				
5	3	D6	0mM glucose 1h	10nM	19:30	20:30		26-35	21-30				
6	4	B6	0mM glucose + FCCP 1h	10nM	19:50	20:50	NA	36-45	31-40				

time table for NCT treatment_MMS_Plate I-BatchB										
	A	B	C	D	E	F	G	H	I	J
1	NCT (plate I)									
2	condition#		condition	drug	media	pyruvate	imaging order	time	imaging	
3										
4	E	C5	MMS 0mM	NT	2.5mM gluc	1mM	2	12:38	12:40	
5	A	C1	MMS 2mM	NT	2.5mM gluc	1mM	4	14:00	15:02	
6	B	C2	MMS 1mM	NT	2.5mM gluc	1mM	6	15:20	16:19	
7	C	C3	MMS 0.5mM	NT	2.5mM gluc	1mM	8	16:40	17:38	position 2 is shifted
8	D	C4	MMS 0.25mM	NT	2.5mM gluc	1mM	10	18:00	18:58	
9										
10	E	D5	MMS 0mM	NT	2.5mM gluc	1mM	1	12:05	12:03	
11	A	D1	MMS 1.5mM	NT	2.5mM gluc	1mM	3	13:25	14:23	
12	B	D2	MMS 1mM	NT	2.5mM gluc	1mM	5	14:40	15:41	
13	C	D3	MMS 0.5mM	NT	2.5mM gluc	1mM	7	16:00	16:58	
14	D	D4	MMS 0.25mM	NT	2.5mM gluc	1mM	9	17:20	18:23	

fileID_treatment_NCT_20260627_MN_110,111						
	A	B	C	D	E	F
1	image	reporter	treatment	drug	time	pyruvate
2	com_110_MN_no-glu_50%_1s_30s_60x_0000	110	0mM MMS	NT	1h	yes
3	com_110_MN_no-glu_50%_1s_30s_60x_0001	110	0mM MMS	NT	1h	yes
4	com_110_MN_no-glu_50%_1s_30s_60x_0002	110	2mM MMS	NT	1h	yes
5	com_110_MN_no-glu_50%_1s_30s_60x_0003	110	2mM MMS	NT	1h	yes
6	com_110_MN_no-glu_50%_1s_30s_60x_0004	110	1mM MMS	NT	1h	yes
7	com_110_MN_no-glu_50%_1s_30s_60x_0005	110	1mM MMS	NT	1h	yes
8	com_110_MN_no-glu_50%_1s_30s_60x_0006	110	0.5mM MMS	NT	1h	yes
9	com_110_MN_no-glu_50%_1s_30s_60x_0007	110	0.5mM MMS	NT	1h	yes
10	com_110_MN_no-glu_50%_1s_30s_60x_0008	110	0.25mM MMS	NT	1h	yes
11	com_110_MN_no-glu_50%_1s_30s_60x_0009	110	0.25mM MMS	NT	1h	yes
12	com_111_MN_no-glu_50%_1s_30s_60x_0000	111	0mM MMS	NT	1h	yes
13	com_111_MN_no-glu_50%_1s_30s_60x_0001	111	0mM MMS	NT	1h	yes
14	com_111_MN_no-glu_50%_1s_30s_60x_0002	111	2mM MMS	NT	1h	yes
15	com_111_MN_no-glu_50%_1s_30s_60x_0003	111	2mM MMS	NT	1h	yes
16	com_111_MN_no-glu_50%_1s_30s_60x_0004	111	1mM MMS	NT	1h	yes
17	com_111_MN_no-glu_50%_1s_30s_60x_0005	111	1mM MMS	NT	1h	yes
18	com_111_MN_no-glu_50%_1s_30s_60x_0006	111	0.5mM MMS	NT	1h	yes
19	com_111_MN_no-glu_50%_1s_30s_60x_0007	111	0.5mM MMS	NT	1h	yes
20	com_111_MN_no-glu_50%_1s_30s_60x_0008	111	0.25mM MMS	NT	1h	yes
21	com_111_MN_no-glu_50%_1s_30s_60x_0009	111	0.25mM MMS	NT	1h	yes

IF on plate I: hypoglycemia

Plate I: IF conditions						
	1	2	3	4	5	6
A						2.5mM gluc
B						0.15mM gluc
C						2.5mM gluc
D						0.15mM gluc

Treat cells at **18:38-43**
Fix cells with 4% PFA for 15min at RT at **19:40**

RNA on plate III: MMS

Well1			
	1	2	3
A	0mM MMS	0.5mM MMS	
B	1mM MMS	2mM MMS	

^

Prepare 2.4mL/condition

Use 1.2mL to wash the well once and add the remaining into the well at around **18:05-18:09**

Collect RNA with 1mL Trizol at **19:05**

2026年6月30日 星期二

IF continue

- 1. 0.2% Triton from 1:21-1:46 (accidentally incubated for a longer time)
- 2. blocking with 10% FBS from 1:46 - 17:55
- 3. Treat 1st Ab

IF list:

RANGAP1 (C-31, Rb)

Importin-b1 (147, M)

antibody												
	A	B	C	D	E	F	G	H	I	J	K	L
1	condition	Ab#	antibody	species	dilution	total volume	Ab volume	secondary Ab	#	dilution	total volume	Ab volume
2	A	C-31	RANGAP1	Rb	500	800	1.6	488	B	1000	800	0.8
3	A	147	Importin-b1	m	1000	800	0.8	594	O	1000	800	0.8

2026年7月1日 星期三

IF continue

- 1. PBS wash*3
- 2. 2ndAb for 1h
- 3. PBS wash*3
- 4. DAPI 5min
- 5. PBS wash*3
- 6. mount

2026年7月3日 星期五

Imaging IF

new Ti2 scope

z stack: 0.5um/step, 10 um in total, 21 steps

NCT_O.C. , 60X, 74X704, 16-bit, old scope						^
	A	B	C	D	E	
1	channel	bit	power (%)	exposure time	note	
2	BF	16		40ms		
3	405	16	15	100 ms		
4	488	16	30	100 ms		
5	561	16	50	600ms	signal is not very bright	

fileID_treatment_IF								/
	A	B	C	D	E	F	G	
1	image	treatment	treat_time	pyruvate	replicate	488	561	
2	A6_488-RANGAP1_561-Impb1_640_405-DAPI_001	2.5mM gluc	1h	1mM pyruvate	A	RanGAP1	Impb1	
3	A6_488-RANGAP1_561-Impb1_640_405-DAPI_002	2.5mM gluc	1h	1mM pyruvate	A	RanGAP1	Impb1	
4	B6_488-RANGAP1_561-Impb1_640_405-DAPI_001	0.15mM gluc	1h	1mM pyruvate	A	RanGAP1	Impb1	
5	C6_488-RANGAP1_561-Impb1_640_405-DAPI_001	2.5mM gluc	1h	1mM pyruvate	B	RanGAP1	Impb1	
6	C6_488-RANGAP1_561-Impb1_640_405-DAPI_002	2.5mM gluc	1h	1mM pyruvate	B	RanGAP1	Impb1	
7	D6_488-RANGAP1_561-Impb1_640_405-DAPI_001	0.15mM gluc	1h	1mM pyruvate	B	RanGAP1	Impb1	