

20260623MN (NCT, ATP, viability, protein)

Project: Yi-Shuan Tseng Holt Lab
Author: Yi Shuan Tseng
Status: In progress
Entry Created On: 15 Jul 2026 15:22:00 UTC
Entry Last Modified: 27 Aug 2026 20:15:00 UTC
Export Generated On: 01 Sep 2026 18:24:14 UTC

2026年6月21日 星期日	(A) T2
2026年6月23日 星期二	(A) T2 -> T3 (B) T2
2026年6月25日 星期四	(A) T3 -> T4 (B) T2 -> T3
2026年6月27日 星期六	(A) T4 -> T5 (B) T3 -> T4
2026年6月29日 星期一	(A) T5 -> T6 (B) T4 -> T5
2026年7月1日 星期三	(A) T6 (B) T5 -> T6
2026年7月3日 星期五	(A) T6 (B) T6
2026年7月6日 星期一	T6 -> T7 Dissociation 1. Collect EBs into 50mL tubes (each plate). Wash the plate with 5mL PBS and added them into the tubes. Let EBs settle with gravity. 2. Remove sup. and add 3mL 0.25% trypsin for each (prepare 8mL 0.25% trypsin with 220U DNaseI) 3. Incubate the tubes in 37C. 4. pipette*10 at 5min, 10min. 5. Stop reaction at ~12 (A), 13 (B) min by adding 8mL FBS into 50mL tube -> pipette with 5mL pipet ~ 10 times 6. add 8mL CTWM -> pipette ~ 10 times 7. filter cells through 40um cell strainer into new tubes. 8. centrifuge at 200xg, 4min

cell info

	A	B	C	D
1		cell	passage	well number
2	A	WT11	P19+2	6
3	B	WT11	P19+3	5

Table1

	A	B
1		volume (mL)
2	trypsin	4
3	HI-FBS/CTWM	8

Table3

	A	B	C	D	E	F	G	H	I	J	K
1		conc. (M/mL)	survival	total volume	tota cell count	plate#	well plate	cell count/well (K)	number of wells	volume	
2	A	3.31	98	5	16.55	24I (A) 92	24	300	8	725	
3	B	3.1	97	5	15.5	24I (B) 92	24	300	8	774	
4						24II (A) 110	24	300	6	544	
5						24II (B) 110	24	300	6	581	
6						24II (A) 111	24	300	4	363	
7						24II (B) 111	24	300	4	387	
8						24II (A) 185	24	300	2	181	
9						24II (B) 185	24	300	2	194	
10						96I (A)	96	100	60	1813	
11						96II (B)	96	100	60	1935	
12						6I (A)	6	1675	2		the rest
13						6I (B)	6	1150	2		the rest

24I

	1	2	3	4	5	6
A	92 A (1uL/well)				A no virus	
B						
C	92 B (1uL/well)				B no virus	
D						

24II

	1	2	3	4	5	6
A	110 A (1uL/well)					
B	111 A (1.6uL/well)				185 A (3uL/well)	
C	110 B (1uL/well)					
D	111 B (1.6uL/well)				185 B (3uL/well)	

96 I and II ^

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

6I ^

	1	2	3
A	A	A	
B	B	B	

2026年7月7日 星期二

Change T7 -> T7

For cells without virus: add more T7

2026年7月9日 星期四

T7 -> T8

2026年7月10日 星期五

1/2 T8 -> 1/2 T8

2026年7月13日 星期一

1/2 T8 -> 1/2 T8

2026年7月15日 星期三

CellTiter Glo

1. Prepare 3000uL media for each and add 230uL in 96 well plate. 1 condition has 6 wells. (6*200+200 = 1400)
2. remove original media with multichannel pipette
3. wash with the media
4. add media and incubate for 1h

Table4			^
	A	B	
1	condition #		
2	1	0mM MMS	
3	2	1mM MMS	
4	3	1.5mM MMS	
5	4	2mM MMS	
6	5	5mM MMS	
7	6	2.5mM gluc	
8	7	0.15mM gluc	
9	8	2.5mM gluc without pyruvate	
10	9	0.15mM gluc without pyruvate	
11	10	0.15mM gluc + FCCP +oligomycin	

Table5			^
	A	B	
1	condition #		
2	1	0mM MMS	
3	2	1mM MMS	
4	3	1.5mM MMS	
5	4	2mM MMS	
6	5	5mM MMS	
7	6	0mM MMS 2h	
8	7	1mM MMS 2h	
9	8	1.5mM MMS 2h	
10	9	2mM MMS 2h	
11	10	5mM MMS 2h	

96I_condition_BatchA													^
	1	2	3	4	5	6	7	8	9	10	11	12	
A													
B	B	1	100m M	150m M			3	100m M	150m M				D
C													
D													
E	C	2					4						E
F													
G													
H	F			G			2.5mM gluc+1mM pyruvate			blank			

96II_condition_BatchB												
	1	2	3	4	5	6	7	8	9	10	11	12
A												
B	B	1	100mM	150mM			4	100mM	150mM			D
C												
D												
E	C	2					5					
F												E
G												
H	F			G			2.5mM gluc+1mM pyruvate			blank		

prepare ATP calibration								/
	A	B	C	D	E	F	G	
1	A	B	C	D	E	F	G	
2	100mM (stock)	1mM	1uM	300nM	100nM	30nM	10nM	
3	dilution factor from the previous one	100	1000	3.333333333	10	3.333333333	10	
4		from A	from B	from C	from C	fromE	fromE	
5	total	1000	1500	1000	1500	1000	1000	
6	added volume	10	1.5	300	150	300	100	
7	media volume	990	1500.0	700	1350	700	900	

prepare media with different conc. of glucose					
	A	B	C	D	E
1	diluted from (mM)	targetconc. (mM)	dilution factor	final volume (uL)	added volume (uL)
2	1100	2.5	440	20000	45.5
3	2.5	0.15	16.6666666667	6000	360.0
4				26000	-> 14mL + 140uL pyruvate



prepare media with different conc. of glucose with pyruvate							
	A	B	C	D	E	F	G
1		diluted from (mM)	targetconc. (uM)	dilution factor	final volume (uL)	added volume (uL)	media volume
2	oligomycin	5	1.5	3333.33333333333	3000	0.9	3000.0
3	FCCP	1	1	1000	3000	3.0	3000.0

drug calculation													
	A	B	C	D	E	F	G	H	I	J	K	L	M
1			stock (mM)	1st dilution factor	total volume for 1st dil.	stock volume	media volume for 1st dil.	conc. after 1st dilution (mM)	working (mM)	dilution factor from 1st dil.	total volume (mL)	volume of 1st dil.	volume of media
2	MMS	6							0		3000	0	3000
3		7	11800	100	100	1	99	118	1	118	3000	25.6	3000
4		8	11800	100	100	1	99	118	1.5	78.6666666667	3000	38.6	3000
5		9	11800	100	100	1	99	118	2	59	3000	51.7	3000
6		10	11800	100	100	1	99	118	5	23.6	3000	132.7	3000
7										total	15000	248.7	

I accidentally used the original MMS for the dilution, so everything is 100X.....I realized it after almost doing the the whole experiment....

I used control experiments for new treatments

wrong concentration ^

	A	B	C	D
1		12:00	prepare media	
2		12:30	prepare media into a 96-well plate	
3	-135	13:09	add media to cells (for 2h treatment)	13:15
4				
5	-75	14:09	add media to cells for 1h treatment	14:15
6	-23	14:55	prepare ATP standard and add samples	
7	-17	15:05	add celltiter glo reagent	15:19
8	-12	15:07	put the plate into cytation 5 and start program to shake	15:21
9	-10	15:09	wait	15:23
10	0	15:19	measure	15;33

redo experiment ^

	A	B	C	D
1		12:00	prepare media	
2	-75	15:10	add media to cells for 1h treatment	15:13
3	-23	15:55	prepare ATP standard and add samples	
4	-17	16:05	add celltiter glo reagent	16:05
5	-12	16:07	put the plate into cytation 5 and start program to shake	16:09
6	-10	16:09	wait/transfer to a new plate	
7	0	16:19	measure	

2026年7月16日 星期四

1/2 T8 -> 1/2 T8 (I forgot to change media yesterday...)

Perceval: BatchA

Plate I Perceval condition							^
	1	2	3	4	5	6	
A	0mM MMS	1mM MMS	1.5m M MMS	0mM MMS + 200nM olaparib			
B	2mM MMS	3mM MMS	5mM MMS	0mM MMS + DMSO			
C	0mM MMS	1mM MMS	1.5m M MMS	2mM MMS + 200nM olaparib			
D	2mM MMS	3mM MMS	5mM MMS	2mM MMS + DMSO			

drug calculation_MMS																		/
	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q	
1														batchA		batchB		
2			stock (mM)	1st dilution factor	total volume for 1st dil.	stock volume	media volume for 1st dil.	conc. after 1st dilution (mM)	working (mM)	dilution factor from 1st dil.	total volume (mL)	volume of 1st dil.	volume of media	treatment time	image time	treatment	image time	
3	MMS	1							0		1000	0	1000	13:28	14:30	21:15	22:19	
4		2	11800	100	100	1	99	118	1	118	1000	8.5	1000	13:55	14:57	21:32	22:35	
5		3	11800	100	100	1	99	118	1.5	78.66666666666667	1000	12.9	1000	14:10	15:13	21:48	22:50	
6		4	11800	100	100	1	99	118	2	59	1000	17.2	1000	14:40	15:40	22:05	23:05	
7		5	11800	100	100	1	99	118	3	39.33333333333333	1000	26.1	1000	15:03	16:02	22:27	23:26	
8		6	11800	100	100	1	99	118	5	23.6	1000	44.2	1000	15:30	16:30	22:43		
9										total	6000	109.0						

fileid_treat_20260709_MN_92						
	A	B	C	D	E	F
1	image	treatment	drug	treat_time	pyruvate	replicate
2	vLH92_MN_60x-1xzoom_0000	0mM MMS	NT	1h	yes	A
3	vLH92_MN_60x-1xzoom_0001	1mM MMS	NT	1h	yes	A
4	vLH92_MN_60x-1xzoom_0002	1.5mM MMS	NT	1h	yes	A
5	vLH92_MN_60x-1xzoom_0003	2mM MMS	NT	1h	yes	A
6	vLH92_MN_60x-1xzoom_0004	3mM MMS	NT	1h	yes	A
7	vLH92_MN_60x-1xzoom_0005	5mM MMS	NT	1h	yes	A
8	vLH92_MN_60x-1xzoom_0006	0mM MMS	NT	1h	yes	B
9	vLH92_MN_60x-1xzoom_0007	1mM MMS	NT	1h	yes	B
10	vLH92_MN_60x-1xzoom_0008	1.5mM MMS	NT	1h	yes	B
11	vLH92_MN_60x-1xzoom_0009	2mM MMS	NT	1h	yes	B
12	vLH92_MN_60x-1xzoom_0010	3mM MMS	NT	1h	yes	B
13	vLH92_MN_60x-1xzoom_0011	5mM MMS	NT	1h	yes	B

non-classical NLS NCT kinetics

time table for NCT treatment_hypoglycemia_Plate II-BatchA								
	A	B	C	D	E	F	G	H
1	NCT (plate I)							
2	condition#		condition	pyruvate	imaging order	time	imaging	
3	110							
4	A	A5	2.5mM gluc	1mM	1	16:05	17:03	XY2 is not usable
5	B	A6	0.15mM gluc	1mM	3	16:35	17:46	
6	A'	C5	0.15mM gluc	1mM	5	18:47	19:46	
7	B'	C6	2.5mM gluc	1mM	7	19:23	20:21	
8	185							
9	A	B5	MMS 0mM	1mM	2	17:25	18:26	
10	B	B6	MMS 1mM	1mM	4	18:05	19:02	
11	A'	D5	MMS 0mM	1mM	6	20:05	21:03	
12	B'	D6	MMS 1mM	1mM	8	20:45	21:43	

fileID_treatment_NCT_20260716_MN_110,185

^

	A	B	C	D	E
1	image	reporter	treatment	time	pyruvate
2	com_110_MN_no-glu_50%_1s_30s_60x_0000	110	2.5mM gluc	1h	yes
3	com_110_MN_no-glu_50%_1s_30s_60x_0001	110	2.5mM gluc	1h	yes
4	com_110_MN_no-glu_50%_1s_30s_60x_0002	110	0.15mM gluc	1h	yes
5	com_110_MN_no-glu_50%_1s_30s_60x_0003	110	0.15mM gluc	1h	yes
6	com_110_MN_no-glu_50%_1s_30s_60x_0004	110	0.15mM gluc	1h	yes
7	com_110_MN_no-glu_50%_1s_30s_60x_0005	110	0.15mM gluc	1h	yes
8	com_110_MN_no-glu_50%_1s_30s_60x_0006	110	2.5mM gluc	1h	yes
9	com_110_MN_no-glu_50%_1s_30s_60x_0007	110	2.5mM gluc	1h	yes
10	com_185_MN_no-glu_50%_1s_30s_60x_0000	185	2.5mM gluc	1h	yes
11	com_185_MN_no-glu_50%_1s_30s_60x_0001	185	2.5mM gluc	1h	yes
12	com_185_MN_no-glu_50%_1s_30s_60x_0002	185	0.15mM gluc	1h	yes
13	com_185_MN_no-glu_50%_1s_30s_60x_0003	185	0.15mM gluc	1h	yes
14	com_185_MN_no-glu_50%_1s_30s_60x_0004	185	0.15mM gluc	1h	yes
15	com_185_MN_no-glu_50%_1s_30s_60x_0005	185	0.15mM gluc	1h	yes
16	com_185_MN_no-glu_50%_1s_30s_60x_0006	185	2.5mM gluc	1h	yes
17	com_185_MN_no-glu_50%_1s_30s_60x_0007	185	2.5mM gluc	1h	yes

drop images_20260716_MN_110,185

^

	A
1	Label
2	com_110_MN_no-glu_50%_1s_30s_60x_0000_XY2

2026年7月17日 星期五

NCT with DMSO and olaparib

time table for NCT treatment_MMS_Plate II-BatchA										
	A	B	C	D	E	F	G	H	I	J
1	NCT (plate I)									
2	condition#		condition	drug	media	pyruvate	imaging order	time	imaging	
3										
4	A	A1	MMS 0mM	DMSO	2.5mM gluc	1mM	1	15:42	16:44	
5	A'	A2	MMS 0mM	olaparib 200nM	2.5mM gluc	1mM	3	17:05	18:04	
6	B	A3	MMS 2mM	DMSO	2.5mM gluc	1mM	5	18:30	19:29	
7	B'	A4	MMS 2mM	olaparib 200nM	2.5mM gluc	1mM	7	19:53	20:55	
8										
9	A	B1	MMS 0mM	DMSO	2.5mM gluc	1mM	2	16:23	17:24	
10	A'	B2	MMS 0mM	olaparib 200nM	2.5mM gluc	1mM	4	17:43	18:43	
11	B	B3	MMS 2mM	DMSO	2.5mM gluc	1mM	6	19:05	20:10	
12	B'	B4	MMS 2mM	olaparib 200nM	2.5mM gluc	1mM	8	20:40	21:42	XY1,2,3,4 have shifted focus...

fileID_treatment_NCT_20260717_MN_110,111

	A	B	C	D	E	F
1	image	reporter	treatment	drug	time	pyruvate
2	com_110_MN_no-glu_50%_1s_30s_60x_0000	110	0mM MMS	DMSO	1h	yes
3	com_110_MN_no-glu_50%_1s_30s_60x_0001	110	0mM MMS	DMSO	1h	yes
4	com_110_MN_no-glu_50%_1s_30s_60x_0002	110	0MM MMS	200nM olaparib	1h	yes
5	com_110_MN_no-glu_50%_1s_30s_60x_0003	110	0mM MMS	200nM olaparib	1h	yes
6	com_110_MN_no-glu_50%_1s_30s_60x_0004	110	2mM MMS	DMSO	1h	yes
7	com_110_MN_no-glu_50%_1s_30s_60x_0005	110	2mM MMS	DMSO	1h	yes
8	com_110_MN_no-glu_50%_1s_30s_60x_0006	110	2mM MMS	200nM olaparib	1h	yes
9	com_110_MN_no-glu_50%_1s_30s_60x_0007	110	2mM MMS	200nM olaparib	1h	yes
10	com_111_MN_no-glu_50%_1s_30s_60x_0000	111	0mM MMS	DMSO	1h	yes
11	com_111_MN_no-glu_50%_1s_30s_60x_0001	111	0mM MMS	DMSO	1h	yes
12	com_111_MN_no-glu_50%_1s_30s_60x_0002	111	0MM MMS	200nM olaparib	1h	yes
13	com_111_MN_no-glu_50%_1s_30s_60x_0003	111	0mM MMS	200nM olaparib	1h	yes
14	com_111_MN_no-glu_50%_1s_30s_60x_0004	111	2mM MMS	DMSO	1h	yes
15	com_111_MN_no-glu_50%_1s_30s_60x_0005	111	2mM MMS	DMSO	1h	yes
16	com_111_MN_no-glu_50%_1s_30s_60x_0006	111	2mM MMS	200nM olaparib	1h	yes
17	com_111_MN_no-glu_50%_1s_30s_60x_0007	111	2mM MMS	200nM olaparib	1h	yes

drop images_20260717_MN_110,111

^

	A
1	Label
2	com_111_MN_no-glu_50%_1s_30s_60x_0006_XY1
3	com_111_MN_no-glu_50%_1s_30s_60x_0006_XY2
4	com_111_MN_no-glu_50%_1s_30s_60x_0006_XY3
5	com_111_MN_no-glu_50%_1s_30s_60x_0006_XY4

time table for Perceval_MMS_Plate I										
	A	B	C	D	E	F	G	H	I	J
1	NCT (plate I)									
2	condition#		condition	drug	media	pyruvate	imaging order	time	imaging	
3										
4	A	A4	MMS 0mM	olaparib 200nM	2.5mM gluc	1mM	1	21:20	22:20	
5	B	B4	MMS 0mM	DMSO	2.5mM gluc	1mM	2	21:45	22:45	
6	A	C4	MMS 2mM	olaparib 200nM	2.5mM gluc	1mM	3	22:05		
7	B	D4	MMS 2mM	DMSO	2.5mM gluc	1mM	4	22:28		

fileid_treat_20260717_MN_92					
	A	B	C	D	E
1	image	treatment	drug	treat_time	pyruvate
2	vLH92_MN_60x-1xzoom_0000	0mM MMS	200nM olaparib	1h	yes
3	vLH92_MN_60x-1xzoom_0001	0mM MMS	DMSO	1h	yes
4	vLH92_MN_60x-1xzoom_0002	2mM MMS	200nM olaparib	1h	yes
5	vLH92_MN_60x-1xzoom_0003	2mM MMS	DMSO	1h	yes

time table for NCT treatment_MMS_Plate II-BatchB										
	A	B	C	D	E	F	G	H	I	J
1	NCT (plate I)									
2	condition#		condition	drug	media	pyruvate	imaging order	time	imaging	
3										
4	A	C1	MMS 0mM	DMSO	2.5mM gluc	1mM	1	13:33	14:33	
5	B	C2	MMS 0mM	olaparib 200nM	2.5mM gluc	1mM	3	14:54	16:57	
6	A'	C3	MMS 2mM	DMSO	2.5mM gluc	1mM	5	16:24	17:25	XY08 had shifted focus
7	B'	C4	MMS 2mM	olaparib 200nM	2.5mM gluc	1mM	7	17:48	18:47	
8										
9	A	D1	MMS 0mM	DMSO	2.5mM gluc	1mM	2	14:12	15:10	cell died a lot. Try to image some alive cells close to the edge, but need to be careful using the data
10	B	D2	MMS 0mM	olaparib 200nM	2.5mM gluc	1mM	4	15:30	16:39	
11	A'	D3	MMS 2mM	DMSO	2.5mM gluc	1mM	6	16:59	18:02	XY01 had shifted focus
12	B'	D4	MMS 2mM	olaparib 200nM	2.5mM gluc	1mM	8	18:23	19:31	delayed 8min

fileID_treatment_NCT_20260718_MN_110,111						
	A	B	C	D	E	F
1	image	reporter	treatment	drug	time	pyruvate
2	com_110_MN_no-glu_50%_1s_30s_60x_0000	110	0mM MMS	DMSO	1h	yes
3	com_110_MN_no-glu_50%_1s_30s_60x_0001	110	0mM MMS	DMSO	1h	yes
4	com_110_MN_no-glu_50%_1s_30s_60x_0002	110	0MM MMS	200nM olaparib	1h	yes
5	com_110_MN_no-glu_50%_1s_30s_60x_0003	110	0mM MMS	200nM olaparib	1h	yes
6	com_110_MN_no-glu_50%_1s_30s_60x_0004	110	2mM MMS	DMSO	1h	yes
7	com_110_MN_no-glu_50%_1s_30s_60x_0005	110	2mM MMS	DMSO	1h	yes
8	com_110_MN_no-glu_50%_1s_30s_60x_0006	110	2mM MMS	200nM olaparib	1h	yes
9	com_110_MN_no-glu_50%_1s_30s_60x_0007	110	2mM MMS	200nM olaparib	1h	yes
10	com_111_MN_no-glu_50%_1s_30s_60x_0000	111	0mM MMS	DMSO	1h	yes
11	com_111_MN_no-glu_50%_1s_30s_60x_0001	111	0mM MMS	DMSO	1h	yes
12	com_111_MN_no-glu_50%_1s_30s_60x_0002	111	0MM MMS	200nM olaparib	1h	yes
13	com_111_MN_no-glu_50%_1s_30s_60x_0003	111	0mM MMS	200nM olaparib	1h	yes
14	com_111_MN_no-glu_50%_1s_30s_60x_0004	111	2mM MMS	DMSO	1h	yes
15	com_111_MN_no-glu_50%_1s_30s_60x_0005	111	2mM MMS	DMSO	1h	yes
16	com_111_MN_no-glu_50%_1s_30s_60x_0006	111	2mM MMS	200nM olaparib	1h	yes
17	com_111_MN_no-glu_50%_1s_30s_60x_0007	111	2mM MMS	200nM olaparib	1h	yes

drop images_20260718_MN_110,111	
	A
1	Label
2	com_110_MN_no-glu_50%_1s_30s_60x_0004_XY8
3	com_111_MN_no-glu_50%_1s_30s_60x_0004_XY1

Viability

Plate 24I viability condition						
	1	2	3	4	5	6
A					0mM MMS+ DMSO	0mM MMS + 200nM olaparib
B					2mM MMS + DMSO	2mM MMS + 200nM olaparib
C					0mM MMS+ DMSO	0mM MMS + 200nM olaparib
D					2mM MMS + DMSO	2mM MMS + 200nM olaparib

Treat cells with SirDNA (1:7000) at 19:08-18
washout & add 1.5uM PI at 20:00
propidium iodide (PI) stock: 1mg/mL (1.5mM)
dilution: 1500uM/1.5uM = 1000
Prepare 1.5uM PI in wash media and add 500uL per well.
Use 20X for imaging

1. Take 9mL 2.5mM gluc media + 1.28uL SirDNA
2. devide media into 4.5uL and add 0.9uL 1/10 DMSO or 200nM olaparib
3. treat cells
4. washout MMS one hour later with media + 1.5uM PI + DMSO/olaparib

**** Images were out of focus...**

2026年7月20日 星期一

Change half media for 6-well plate

2026年7月21日 星期二

Protein extraction

I want to check is the level of SUMOylated RANGAP1 is differernt.

Well1			
	1	2	3
A	2.5mM gluc	0.15mM gluc	
B	2.5mM gluc	0.15mM gluc	

1. Wash cells with 1mL media and treat them at 16:22
2. Precool centrifuger and eppendorf

3. Prepare RIPA buffer with 2%SDS (250uL RIPA + 25uL 20%SDS)
4. Incubate cells in 37C incubator
5. At 17:22, placed the whole plate on ice. And process the following steps one well by one well.
 - a. remove media carefully with P1000
 - b. took 500uL cold PBS and flush the cells down quickly and strongly to detach the cells
 - c. Collect cells suspension into the pre-cooled eppendorf
 - d. After process every cells, I rinse the well again with 500uL cold PBS
 - e. ** I checked the plate under microscope and most cells were detached. Some cells were still there, but I saw many non-MNs are still attached on the well!!! It can be a good way to get rid of non-MNs.
6. centrifuge the cells with 500xg, 2min
7. Add 2.5 µl PMSF solution, 2.5 µl sodium orthovanadate solution and 2.5 µl protease inhibitor cocktail solution (1:100) freshly
8. remove supernatant with P1000 and leave a bit PBS to avoid aspirate out the cell pellet
9. estimate the amount of the left liquid and add approximately same amount of RIPA. I added 40uL this time to not dilute protein too much. Pipette the cells a bit.
10. Boil the lysate at 95C for 5-10min
11. Store the protein in -80C freezer