

20260831EB_MN

Project: Yi-Shuan Tseng Holt Lab

Author: Yi Shuan Tseng

Status: In progress

Entry Created On: 31 Aug 2026 20:12:52 UTC

Entry Last Modified: 31 Aug 2026 20:13:18 UTC

Export Generated On: 01 Sep 2026 18:26:10 UTC

2026年8月31日 星期一

split iPSCs into a full 6-well plate on the Friday before differentiation
differentiate 5 of 6-wells with 10mL T2 in cell-repellent plates
keep some iPSCs

2026年9月2日 星期三

T2 -> T3

2026年9月4日 星期五

T3 -> T4 (~14 mL for over-weekend)

2026年9月7日 星期一

T4 -> T5

2026年9月9日 星期三

T5 -> T6

2026年9月11日 星期五

T6 -> T6 (~14mL for over-weekend)

2026年9月14日 星期一

T6 -> T7

Dissociation

1. Collect EBs into 50mL tubes (each plate). Wash the plate with 5mL PBS and added them into the tubes. Let EBs settle with gravity.
2. Remove sup. and add 3mL 0.25% trypsin for each (prepare 8mL 0.25% trypsin with 220U DNaseI)
3. Incubate the tubes in 37C.
4. pipette*10 at 5min, 10min.
5. Stop reaction by adding 8mL FBS into 50mL tube -> pipette with 5mL pipet ~ 10 times
6. add 8mL CTWM -> pipette ~ 10 times
7. filter cells through 40um cell strainer into new tubes.
8. centrifuge at 200xg, 4min

cell info ^

	A	B	C
1	cell	passage	well#
2	WT11		

Table1 ^

	A	B
1		volume (mL)
2	trypsin	3
3	HI-FBS/CTWM	6

Table3 ✓

	A	B	C	D	E	F	G	H	I
1	conc. (M/mL)	survival	total volume	tota cell count	plate#	well plate	cell count/well (K)	number of wells	volume
2	4.26	99	7	29.82	A	96	60	36	507
3					B	24	250	16	939
4					C,D	6	1972	12	the rest

2026年9月15日 星期二

Change T7 -> T7 if added virus

2026年9月17日 星期四

T7 -> T8

2026年9月18日 星期五

1/2 T8 -> 1/2 T8

2026年9月21日 星期一

1/2 T8 -> 1/2 T8

2026年9月23日 星期三

1/2 T8 -> 1/2 T8