

RanGAP1 SUMOylation

Project: Yi-Shuan Tseng Holt Lab

Author: Yi Shuan Tseng

Status: In progress

Entry Created On: 28 Aug 2026 22:02:03 UTC

Entry Last Modified: 01 Sep 2026 17:57:22 UTC

Export Generated On: 01 Sep 2026 18:23:04 UTC

2026年8月28日 星期五

Goal: I want to check if hypoglycemia changes SUMOylation level of RanGAP1

Cells: 20260623MN (4 samples) and 20260629MN (10 samples)

MW:

- unmodified RanGAP1: ~65kDa
- SUMOylated RanGAP1: ~82kDa

2026年9月1日 星期二

BCA assay

image.png

Table 1 BSA standard dilution scheme for the microplate procedure and the test tube procedure using the standard protocol (working range = 20-2,000 µg/mL)

Vial	Volume of diluent	Volume and source of BSA	Final BSA concentration
A	0	300 µL of stock	2,000 µg/mL
B	125 µL	375 µL of stock	1,500 µg/mL
C	325 µL	325 µL of stock	1,000 µg/mL
D	175 µL	175 µL of vial B dilution	750 µg/mL
E	325 µL	325 µL of vial C dilution	500 µg/mL
F	325 µL	325 µL of vial E dilution	250 µg/mL
G	325 µL	325 µL of vial F dilution	125 µg/mL
H	400 µL	100 µL of vial G dilution	25 µg/mL
I	400 µL	0	0 = Blank

Well1													^
	1	2	3	4	5	6	7	8	9	10	11	12	
A	A			1			9						
B	B			2			10						
C	C			3			11						
D	D			4			12						
E	E			5			13						
F	F			6			14						
G	G			7			blank						
H	H			8									

Working reagent

(9 standard + 14 unknown)*3*200=13800

14mL A + 280uL B