

3D neruon organoid with microfluidic device

Project: Yi-Shuan Tseng Holt Lab
Author: Yi Shuan Tseng
Status: In progress
Entry Created On: 15 Jun 2026 23:13:20 UTC
Entry Last Modified: 15 Jun 2026 23:19:01 UTC
Export Generated On: 01 Sep 2026 18:34:29 UTC

2026年6月8日 星期一	Yoshiho attached a microfluidic device on a coverslip with ethanol and let it dry out
2026年6月9日 星期二	Coat a half of device with matrigel in stemflex and another half with PLO for 30 min. Seed one EB into each chamber (20260525EB_MN)
2026年6月10日 星期三	Chambers coated with PLO were contaminated, but those with matrigel were fine. (at the end I found out the PBS was contaminated) Add 1M NaOH into the contaminated chambers for a few minutes and aspirate the solution out to save the remaining wells. Cells with matrigel grew a bit neurites
2026年6月11日 星期四	Cells with matrigel survived well
2026年6月15日 星期一	